MRIMath

Generated by Doxygen 1.8.15

Contents

1	Nam	espace Index	1
	1.1	Namespace List	1
2	Clas	s Index	3
	2.1	Class List	3
3	Nam	espace Documentation	5
	3.1	DataHandler Namespace Reference	5
		3.1.1 Detailed Description	5
	3.2	EmailHandler Namespace Reference	5
		3.2.1 Detailed Description	5
	3.3	GeneralNetwork Namespace Reference	5
		3.3.1 Detailed Description	5
	3.4	HardwareHandler Namespace Reference	6
		3.4.1 Detailed Description	6
	3.5	LoadAndTestModel Namespace Reference	6
		3.5.1 Detailed Description	6
	3.6	TimerModule Namespace Reference	6
		3.6.1 Detailed Description	7
	3.7	TrainModels Namespace Reference	7
		3.7.1 Detailed Description	7

ii CONTENTS

4	Clas	s Docu	mentation	9
	4.1	DataHa	andler.DataHandler Class Reference	9
		4.1.1	Constructor & Destructor Documentation	10
			4.1.1.1init()	10
		4.1.2	Member Function Documentation	10
			4.1.2.1 derivePatches()	10
			4.1.2.2 derivePatchFromSegments()	11
			4.1.2.3 deriveRandomPatch()	11
			4.1.2.4 getData()	12
			4.1.2.5 getDirectoryFromIndex()	12
			4.1.2.6 getImage()	12
			4.1.2.7 loadDataParallel()	13
			4.1.2.8 loadDataSequential()	13
			4.1.2.9 loadIndividualPatient()	13
	4.2	Email	landler.EmailHandler Class Reference	14
		4.2.1	Constructor & Destructor Documentation	14
			4.2.1.1init()	15
		4.2.2	Member Function Documentation	15
			4.2.2.1 attachFile()	15
			4.2.2.2 prepareMessage()	15
			4.2.2.3 sendMessage()	15
		4.2.3	Member Data Documentation	16
			4.2.3.1 addressBook	16
	4.3	Hardwa	areHandler.HardwareHandler Class Reference	16
		4.3.1	Constructor & Destructor Documentation	17
			4.3.1.1init()	17
		4.3.2	Member Function Documentation	17
			4.3.2.1 createThreadPool()	17
			4.3.2.2 getAvailableGPUs()	17
			4.3.2.3 getNumberOfCores()	18
	4.4	TimerN	Module.TimerModule Class Reference	18
		4.4.1	Member Function Documentation	18
			4.4.1.1 getElapsedTime()	18
Inc	dex			19

Chapter 1

Namespace Index

1.1 Namespace List

Here is a list of all documented namespaces with brief descriptions:

DataHandler	5
EmailHandler	5
GeneralNetwork	5
HardwareHandler	6
LoadAndTestModel	
TimerModule	
TrainModels	7

2 Namespace Index

Chapter 2

Class Index

2.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

DataHandler.DataHandler	9
EmailHandler.EmailHandler	14
HardwareHandler.HardwareHandler	16
TimerModule.TimerModule	18

4 Class Index

Chapter 3

Namespace Documentation

3.1 DataHandler Namespace Reference

Classes

class DataHandler

3.1.1 Detailed Description

Class designed to do all data handling and manipulation, ranging from dataloading to network preprocessing. As time goes on, some of this may be refactored, and some of this functionality is contingent on the data being in a certain structure.

@author Daniel Enrico Cahall

3.2 EmailHandler Namespace Reference

Classes

· class EmailHandler

3.2.1 Detailed Description

Class designed to handle constructing and sending emails, usually in the context of notifying one or more recipients when a process has finished. Currently, the EmailHandler has the capability to send an email to one or more people as long as they are identified in the Address Book, and attach one or more files to the email

@author Daniel Enrico Cahall

3.3 GeneralNetwork Namespace Reference

3.3.1 Detailed Description

Created on Jan 12, 2018

@author: daniel

3.4 HardwareHandler Namespace Reference

Classes

· class HardwareHandler

3.4.1 Detailed Description

Class designed to handle all hardware related tasks, such as creating a threadpool or getting the number of co I anticipate that this class will grow over time, but for the time being it handles all necessary hardware tas

3.5 LoadAndTestModel Namespace Reference

Variables

- string **data_dir** = '/media/daniel/ExtraDrive1/Patient_Data_Images';
- model = load model('/home/daniel/eclipse-workspace/MRIMath/Models/2018-01-25 20 55/model.h5');
- dataHandler = DataHandler()
- list segments = []
- patient_directory = os.fsencode(dataHandler.getDirectoryFromIndex(i, data_dir))
- string **orig_data** = patient_directory + b'/Original_Img_Data'
- string **seg_data** = patient_directory + b'/Segmented_Img_Data'
- **img** = dataHandler.getImage(orig_data+b'/'+file)
- patches = dataHandler.derivePatches(img, 1)
- x
- у
- **patch** = patch.reshape(1,25, 25, 1)
- **pred** = model.predict(patch)
- **label** = np.argmax(pred)
- **fig** = plt.figure()
- **a** = fig.add_subplot(1,9,1)
- int **ind** = 2;

3.5.1 Detailed Description

```
Created on Nov 28, 2017
@author: daniel
```

3.6 TimerModule Namespace Reference

Classes

· class TimerModule

3.6.1 Detailed Description

Class designed to keep track of time and performance. Pretty simple and small, although more capability can be if necessary.

@author Daniel Enrico Cahall

3.7 TrainModels Namespace Reference

Functions

• def precision (y_true, y_pred)

Variables

- **now** = datetime.now()
- date_string = now.strftime('%Y-%m-%d %H %M')
- dataHandler = DataHandler()
- emailHandler = EmailHandler()
- hardwareHandler = HardwareHandler()
- timer = TimerModule()
- tuple input_img = (dataHandler.n, dataHandler.n, 1)
- model = Sequential()
- string data_dir = '/coe_data/MRIMath/MS_Research/Patient_Data_Images'
- · training
- training_labels
- · testing
- · testing_labels
- string model directory = "/coe data/MRIMath/MS Research/MRIMath/Models/" + date string
- int **num_epochs** = 15
- int batchSize = 64
- float **Irate** = 0.1
- float momentum = 0.9
- **sgd** = SGD(Ir=Irate, momentum=momentum, nesterov=True)
- string model_info_filename = 'model_info.txt'
- model_info_file = open(model_directory + '/' + model_info_filename, "w")
- string log_info_filename = 'model_loss_log.csv'
- log_info = open(model_directory + '/' + log_info_filename, "w")
- csv_logger = CSVLogger(model_directory + '/' + log_info_filename, append=True, separator=',')
- reduce Ir = ReduceLROnPlateau(monitor='val loss', factor=0.2, patience=3, min lr=0.001)
- **G** = hardwareHandler.getAvailableGPUs()
- parallel_model = multi_gpu_model(model, G)
- · optimizer
- loss
- metrics
- · epochs
- · batch_size
- shuffle
- · validation data
- · callbacks
- string message = "Finished training network at " + str(datetime.now()) + '\n\n'
- · print fn

3.7.1 Detailed Description

Created on Jan 9, 2018 @author: daniel

Chapter 4

Class Documentation

4.1 DataHandler.DataHandler Class Reference

Public Member Functions

• def __init__ (self, tolerance=0.25, numPatches=10, n=25)

The constructor for the datahandler class.

def getImage (self, path)

Reads an image from a filepath.

• def loadDataSequential (self, data_directory, start, finish)

Loads patient images and segments sequentially, assuming you want to through a range of numbered patients.

• def loadDataParallel (self, data_directory, start, finish)

Loads patient images and segments in parallel, assuming you want to through a range of numbered patients.

def loadIndividualPatient (self, index, data_directory)

Derives and labels patches from an individual patient.

• def getDirectoryFromIndex (self, index, data_directory)

Constructs the patient directory string based on the index (based on current labeling scheme)

def preprocessForNetwork (self)

Preprocesses the data for the network by converting the list of patches and labels to a numpy array, and normalizing the patches.

• def deriveRandomPatch (self, patient_directory, img, file)

Derives random patches from an image.

• def derivePatches (self, img, stepSize)

Derives an individual patch from an image.

def derivePatchFromSegments (self, patient_dir, x, y, img_num)

Derives and labels the patches in the segment imiage.

def getData (self)

Acquires the data from the DataHandler after all data has been loaded and processed.

def clearVectors (self)

Clears the data and labels.

Public Attributes

- X
- labels

Static Public Attributes

```
• lock = threading.Lock()
```

- manager = Manager()
- list **X** = []
- list **labels** = []
- int **W** = 240
- int **H** = 240
- hardwareHandler = HardwareHandler()
- tolerance = None
- numPatches = None
- **n** = None
- list patches = []
- **window** = img[y:y + self.n, x:x + self.n]

4.1.1 Constructor & Destructor Documentation

The constructor for the datahandler class.

Parameters

tolerance	the percentage of pixels in a patch that can be background (default 0.25)
numPatches	the number of patches to extract per image (default 10)
n	the dimensions of the patch to be taken from the image (default 25)

4.1.2 Member Function Documentation

4.1.2.1 derivePatches()

```
def DataHandler.DataHandler.derivePatches ( self, \\ img, \\ stepSize )
```

Derives an individual patch from an image.

Parameters

img	the image to derive patches from
stepSize	the amount the sliding window shifts per iteration
file	the patient image number (e.g. img_1)

Returns

patches a list of all patches in the image

4.1.2.2 derivePatchFromSegments()

Derives and labels the patches in the segment imiage.

Parameters

patient_dir	the specific patient directory (e.g. Patient_001_Data)
X	the starting point for columns (x-direction) for the patch
У	the starting point for rows (y-direction) for the patch
img_num	image number (e.g. img_1)

Returns

a boolean flag which states if a label for the segment was succesfully found

4.1.2.3 deriveRandomPatch()

Derives random patches from an image.

Parameters

patient_directoy	the directory where the specific patient data is located (e.g. Patient_001_Data)
img	the image to derive patches from
Gefikeated by Doxygen	the patient image number (e.g. img_1)

4.1.2.4 getData()

```
\begin{tabular}{ll} \tt def DataHandler.DataHandler.getData ( \\ & self ) \end{tabular}
```

Acquires the data from the DataHandler after all data has been loaded and processed.

Parameters

```
img_num image number (e.g. img_1)
```

Returns

the data and the labels for the loaded and processed data

4.1.2.5 getDirectoryFromIndex()

```
def DataHandler.DataHandler.getDirectoryFromIndex ( self, \\ index, \\ data\_directory \; )
```

Constructs the patient directory string based on the index (based on current labeling scheme)

Parameters

index	the index of the patient that you need the specific directory for
data_directory	Directory where all patient data is located

4.1.2.6 getImage()

Reads an image from a filepath.

Parameters

path	the path to an image file	
patit	ino patri to ari imago mo	

Returns

the image from the filepath (if one existed) as a numpy array

4.1.2.7 loadDataParallel()

Loads patient images and segments in parallel, assuming you want to through a range of numbered patients.

Parameters

data_directory	the directory where all patient data is located
start	the patient number to start with (inclusive)
finish	the patient number to stop at (exclusive)

4.1.2.8 loadDataSequential()

Loads patient images and segments sequentially, assuming you want to through a range of numbered patients.

Parameters

data_directory	the directory where all patient data is located	
start	the patient number to start with (inclusive)	
finish	the patient number to stop at (exclusive)	

4.1.2.9 loadIndividualPatient()

Derives and labels patches from an individual patient.

Parameters

data_directory	the directory where all patient data is located	
index	index of the patient in the numbered directory	

The documentation for this class was generated from the following file:

· DataHandler.py

4.2 EmailHandler.EmailHandler Class Reference

Public Member Functions

def __init__ (self)

The constructor for the emailHandler class.

• def prepareMessage (self, subject, body)

Prepares the message to be sent by setting up the subject and body.

def connectToServer (self)

Connects to the gmail server and logs in using the mrimathnotifier gmail address.

• def sendMessage (self, recipients)

Sends the email to all desired recipients as long as they are within the address book.

• def finish (self)

Clears the body of the email, the attached files, and the list of recipients, and disconnects from the gmail server.

• def attachFile (self, file, filename)

Attaches a file to the email.

Public Attributes

- msg
- body
- server

Static Public Attributes

- string addr = "mrimathnotifier@gmail.com"
- string **password** = "mrimathpw"
- · dictionary addressBook
- string body = ""

4.2.1 Constructor & Destructor Documentation

The constructor for the emailHandler class.

This creates the message and sets the from address to the mrimathnotifier@gmail.com

4.2.2 Member Function Documentation

4.2.2.1 attachFile()

Attaches a file to the email.

Parameters

file	the file to attach to the email
filename	the name of the file to attach (may not be necessary actually)

4.2.2.2 prepareMessage()

Prepares the message to be sent by setting up the subject and body.

Parameters

subject	Subject line of the email to be sent	
body	Body of the email to be sent	

4.2.2.3 sendMessage()

```
{\tt def\ Email Handler. Email Handler. send Message\ (}
```

```
self,
recipients )
```

Sends the email to all desired recipients as long as they are within the address book.

Parameters

recipients | a list of the names of desired recipients which have their emails linked in the address book map

4.2.3 Member Data Documentation

4.2.3.1 addressBook

```
dictionary EmailHandler.EmailHandler.addressBook [static]
```

Initial value:

The documentation for this class was generated from the following file:

· EmailHandler.py

4.3 HardwareHandler.HardwareHandler Class Reference

Public Member Functions

```
• def init (self)
```

The constructor for the hardwarehandler class.

• def getAvailableGPUs (self)

Acquires the number of GPUs available to use.

• def getNumberOfCores (self)

Acquires the number of CPU cores to use.

def createThreadPool (self, threadCount=None)

Creates a threadpool, where the number of threads is the number of available cores by default.

Public Attributes

numThreads

Static Public Attributes

• pool = None

4.3.1 Constructor & Destructor Documentation

The constructor for the hardwarehandler class.

Sets the default number of threads to the number of available cores

4.3.2 Member Function Documentation

4.3.2.1 createThreadPool()

Creates a threadpool, where the number of threads is the number of available cores by default.

Parameters

threadCount	the number of threads to use (default is number of cores)
-------------	---

Returns

pool the threadpool which was created

4.3.2.2 getAvailableGPUs()

```
\label{eq:continuous} \mbox{ def HardwareHandler.HardwareHandler.getAvailableGPUs (} \\ self \mbox{ )}
```

Acquires the number of GPUs available to use.

Returns

the available number of GPUs on the device (int)

4.3.2.3 getNumberOfCores()

```
\label{eq:continuous} \mbox{def HardwareHandler.HardwareHandler.getNumberOfCores (} \\ self \mbox{)}
```

Acquires the number of CPU cores to use.

Returns

the number of cores on the device (int)

The documentation for this class was generated from the following file:

HardwareHandler.py

4.4 TimerModule.TimerModule Class Reference

Public Member Functions

def startTimer (self)

Starts the timer.

def stopTimer (self)

Stops the timer.

def getElapsedTime (self)

Computes the amount of time elapsed based on when the timer started and stopped.

Public Attributes

- · start_time
- · stop_time

Static Public Attributes

- float start_time = 0.0
- float stop_time = 0.0

4.4.1 Member Function Documentation

4.4.1.1 getElapsedTime()

```
\label{eq:continuous} \mbox{def TimerModule.TimerModule.getElapsedTime (} \\ self \mbox{)}
```

Computes the amount of time elapsed based on when the timer started and stopped.

Returns

the amount of time between the time being started and stopped

The documentation for this class was generated from the following file:

· TimerModule.py

Index

init DataHandler::DataHandler, 10 EmailHandler::EmailHandler, 14 HardwareHandler::HardwareHandler, 17	getImage DataHandler::DataHandler, 12 getNumberOfCores HardwareHandler::HardwareHandler, 17
addressBook EmailHandler::EmailHandler, 16 attachFile EmailHandler::EmailHandler, 15 createThreadPool HardwareHandler::HardwareHandler, 17 DataHandler, 5 DataHandler.DataHandler, 9 DataHandler::DataHandlerinit, 10	HardwareHandler, 6 HardwareHandler.HardwareHandler, 16 HardwareHandler::HardwareHandlerinit, 17 createThreadPool, 17 getAvailableGPUs, 17 getNumberOfCores, 17 LoadAndTestModel, 6 loadDataParallel DataHandler::DataHandler, 13 loadDataSequential
derivePatchFromSegments, 11 derivePatches, 10 deriveRandomPatch, 11 getData, 12 getDirectoryFromIndex, 12 getImage, 12 loadDataParallel, 13 loadDataSequential, 13 loadIndividualPatient, 13	DataHandler::DataHandler, 13 loadIndividualPatient DataHandler::DataHandler, 13 prepareMessage EmailHandler::EmailHandler, 15 sendMessage EmailHandler::EmailHandler, 15
derivePatchFromSegments DataHandler::DataHandler, 11 derivePatches DataHandler::DataHandler, 10 deriveRandomPatch DataHandler::DataHandler, 11	TimerModule, 6 TimerModule.TimerModule, 18 TimerModule::TimerModule getElapsedTime, 18 TrainModels, 7
EmailHandler, 5 EmailHandler.EmailHandler, 14 EmailHandler::EmailHandlerinit, 14 addressBook, 16 attachFile, 15 prepareMessage, 15 sendMessage, 15	
GeneralNetwork, 5 getAvailableGPUs HardwareHandler::HardwareHandler, 17 getData DataHandler::DataHandler, 12 getDirectoryFromIndex DataHandler::DataHandler, 12 getElapsedTime TimerModule::TimerModule, 18	